



BIOLABO
www.biolabo.fr

MANUFACTURER:
BIOLABO SAS,
Les Hautes Rives
02160, Maizy, France

UREA U.V Kinetic Method

Reagent for quantitative determination of urea
in human serum and plasma or urines.

REF	92032	R1	7 x 30 mL	R2	7 x 30 mL	R3	1 x 10 mL
REF	92132	R1	10 x 100 mL	R2	10 x 100 mL	R3	1 x 10 mL



Made In France

TECHNICAL SUPPORT AND ORDERS

Tel: (33) 03 23 25 15 50

support@biolabo.fr

Latest revision: www.biolabo.fr

I: corresponds to significant modifications

I INTENDED USE

This reagent is designated for professional use in laboratory (automated method).


It allows the quantitative determination of urea (UREA) in human serum and plasma or urines.

I GENERALITIES (1) (6)

More than 90% of urea is excreted through the kidneys in urines. Measurement of the plasma or serum urea concentration is widely regarded as a test of renal function. However, a number of non-renal factors also influence the circulating urea concentration: Urea increased level occurs when proteins catabolism is accelerated, burns, stress, myocardial infarction... Urea is decreased in acute liver destruction and is accompanied with increased ammonium level. Urea level is generally studied in conjunction with creatinine level (urea/creatinine ratio) to refine the diagnosis of post-renal or pre-renal azotemia.

PRINCIPLE (4) (5)

Enzymatic method based on Talke and Schubert reaction, simplified by Tiffany and al. who demonstrated that urea concentration is proportional to absorbance change at 340 nm over a fixed time interval. Reaction scheme is as follows:



REAGENTS

R1 UREA BUF ENZ Buffer Enzymes

Tris pH 7.9 ± 0.1 at 30°C 80 mmol/L

Oxoglutarate 5 mmol/L

Preservative

R2 UREA COENZ Coenzyme

NADH ≥ 0,2 mmol/L

Urease 20000 IU/L

GLDH ≥ 1200 IU/L

R3 UREA Standard

Urea 40 mg/dL (6.66 mmol/L)

According to 1272/2008/EC Regulation, these reagents are not classified as dangerous.

SAFETY CAUTIONS

- Refer to current Material Safety Data Sheet available on request or on www.biolabo.fr
- Verify the integrity of the contents before use.
- Waste disposal: Respect legislation in force in the country.
- All specimens or reagents of biological origin should be handled as potentially infectious. Respect legislation in force in the country.

I Any serious incident that has occurred in connection with the device is notified to the manufacturer and the competent authority of the Member State in which the user and/or patient is based.

REAGENTS PREPARATION

Use a non-sharp instrument to remove aluminum cap.

Add promptly the content of vial R2 into vial R1.

Mix gently until complete dissolution.

Vial R3: Ready to use

STABILITY AND STORAGE

Stored away from light, well capped in the original vial at 2-8°C, reagents are stable when stored and used as described in the insert:

Unopened:

- Until expiry date stated on the label of the kit

Once opened:

- Working reagent is stable for 1 month when free from contamination
- Discard reagent if cloudy or if absorbance of working reagent measured at 340 nm is < 1.100.

SPECIMEN COLLECTION AND HANDLING (2)

Unhemolysed serum or heparinised plasma.

Avoid fluoride or ammonium as anticoagulants which interfere with the assay.

- Stable for 24 h at room temperature
- several days at 2-8°C
- at least 2-3 months frozen

24h Urines:

- Stable for 4 days at 2-8°C
- Add antibacterial agent as Thymol to improve the stability
- Dilute (1+19) with demineralised water before assay

LIMITS (3)

For a more comprehensive review of factors affecting this assay refer to the publication of Young D.S.

MATERIAL REQUIRED BUT NOT PROVIDED

1. Medical analysis laboratory equipment.
2. Spectrophotometer or Biochemistry Clinical Analyzer

Manufacturer	Expiry date	In vitro diagnostic	Storage temperature	Dematerialized water	Biological risk
Product Reference	See Insert	Batch number	Store away from light	Sufficient for	Dilute with

QUALITY CONTROL

- **REF** 95010 EXATROL-N Level I
- **REF** 95011 EXATROL-P Level II
- **REF** 95012 Urinary Controls
- External quality control program

It is recommended to control in the following cases:

- At least once a run
- At least once within 24 hours
- When changing vial of reagent
- After maintenance operations on the instrument

If control is out of range, apply following actions:

1. Prepare a fresh control serum and repeat the test
 2. If control is still out of range, use a new vial of fresh calibrator
 3. If control is still out of range, use a new vial of reagent and reassay
- If control is still out of range, please contact BIOLABO technical support or your local Agent.

EXPECTED VALUES (2)

Serum or plasma	(mg/dL)	[mmol/L]
In cord	45-86	[7.5-14.3]
Premature	6-54	[1.1-8.9]
< 1 year	9-41	[1.4-6.8]
Children	11-39	[1.8-6.4]
18-60 years	13-43	[2.1-7.1]
60-90 years	17-49	[2.9-8.2]
> 90 years	21-66	[3.6-11.1]
Urines	26-43 g/24 h	[0.43-0.71 mol/24 h]

Each laboratory should establish its own normal ranges for the population that it serves.

PERFORMANCES

With Procedure n°1:

On Cobas Mira, at 37°C, 340 nm:

Linearity range: from 20 to 250 mg/dL

Precision:

Within-run N = 20	Low level	Normal level	High level	Between run N = 20	Low level	Normal level	High level
Mean (mg/dL)	23	60	146	Mean (mg/dL)	16	39	131
S.D. mg/dL	0.48	1.22	2.25	S.D. mg/dL	0.46	0.90	2.1
C.V. %	2.07	2.04	1.54	C.V. %	2.9	2.3	1.6

Comparison studies with commercially available reagent:

Realised on human specimens (n=93) between 20 and 250 mg/dL
 $y = 0.9961x + 0.16$ $r = 0.9970$

Interferences:

Ascorbic acid	No interference up to 2500 mg/dL
Total bilirubin	No interference up to 583 µmol/L
Haemoglobin	No interference up to 248 µmol/L
Glucose	No interference up to 1110 mg/dL
Turbidity	Positive interference from 0.320 OD

Other substances may interfere (see § Limits)

Calibration stability :

Make a new calibration when changing reagent batch, if quality control results are found out of the established range and after maintenance operations.

CALIBRATION (7)

- **REF** 95015 Multicalibrator traceable to SRM 909

Or

- Standard (vial R3) for manual procedure only

The calibration frequency depends on proper instrument functions and on the preservation of reagent

PROCEDURE

Manual method:

Let stand reagents and specimens at room temperature.

Procedure n°1

Pipette into thermostated cuvette (37°C):	Standard	Assay
Reagent	1 mL	1 mL
Standard	5 µL	
Specimen (1)		5 µL

Mix. Start a timer. After 30 seconds, record absorbance A1 at 340 nm and then absorbance A2 after 90 seconds.

Procedure n°2

Pipette into thermostated cuvette (37°C):	Standard	Assay
Reagent	1 mL	1 mL
Standard	10 µL	
Specimen (1)		10 µL

Mix. Start a timer. After 30 seconds, record absorbance A1 at 340 nm and then absorbance A2 after 90 seconds.

1. Serum, plasma or urines diluted (1+19) in demineralised water.
2. Performances with manual procedure should be validated by user.
3. KENZA applications and other proposal are available on request.

CALCULATION

Manual Procedure:

- Serum and plasma:

$$\text{Result} = \frac{\text{Abs (A1 - A2) Assay}}{\text{Abs (A1 - A2) Standard}} \times \text{Standard concentration}$$

To calculate blood urea nitrogen (BUN): multiply the value of urea (mg/dL) by 0.467.

- Urines diluted (1+19):

Multiply the result by 20 (dilution factor).

Automatic Biochemistry analyzer:

The analyzer provides directly final result.

For more details about calibration and calculation of results, refer to User's manual and specific application.

REFERENCES

- (1) TIETZ N.W. *Text book of clinical chemistry*, 3rd Ed. C.A. Burtis, E.R. Ashwood, W.B. Saunders (1999) p. 1239-1241.
- (2) *Clinical Guide to Laboratory Test*, 4th Ed., N.W. TIETZ (2006) p. 1096-1099.
- (3) YOUNG D.S., *Effect of Drugs on Clinical laboratory Tests*, 4th Ed. (1990) p. 3-599 to 3-609
- (4) Talke H., Schubert G. E., *Klin. Wochschr.*, 19, (1965), 43, p.174
- (5) Tiffany T. O., and al., *Clin. Chem.*, 18, (1972) p.829-840
- (6) Bernard S. *Bioch. clin. Diagnostics médicaux chirurgicaux 2^{ème} éd.* p.143 144. Ed. Maloine PARIS (1989).
- (7) SRM: Standard Reference Material ®